

journey SoftSecure

ROTATING TRANSFER BENCH

PRODUCT MANUAL



ROTATING TRANSFER BENCH

 journey
HEALTH & LIFESTYLE®

Welcome to Journey Health & Lifestyle!

Thank you for purchasing your SoftSecure Rotating Transfer Bench. You've invested in a high-quality, premium materials rotating transfer bench that will facilitate comfort, safety and functionality while you get in and out of your tub.

Safety instructions

- Read these instructions for use carefully.
- All components should be checked for damage and a secure fit prior to use.
- Approved for indoor use only.
- Clothing or body parts can be pinched when sitting down or standing up, as well as when raising the armrests and locking them into place.
- Prevent unauthorized use, for example, by children.
- Take note of the maximum weight permitted.

Parts of the product



1. Flip Up Armrest
2. Seat
3. Height adjustable leg
4. Rubber tip
5. Backrest
6. Rotating System
7. Sliding Lock
8. Side frame

Assembly

Step 1

Attach the side frames by sliding the tubing into the main section receiver and securing them with the knobs provided.


Holes on main frame need to face outward in order for sliding mechanism to function.

Hole cannot be covered. It will prevent product from sliding.



Step 2

Attach backrest to the frame by depressing push buttons and sliding the backrest into the receivers on the rear sides of the seat.



Step 3

Adjust legs to your height requirement.



Step 4

To flip up armrest, press button on side of seat.

III. INSTRUCTIONS FOR USE

Step 5

- (1) To rotate seat, pull out knob and turn seat until it locks into position.
- (2) Unlock the seat by pulling out the knob in the middle of the rotating mechanism.
- (3) When using transfer bench, ensure seat is locked into place.



Step 6

To slide the seat side-to-side, pull-out red knob.



Product Specifications	
Product Dimension	(33"-34") x (19"-20") x (32"-36")
Weight Capacity	300 lbs
Packaging Dimensions	27" x 18" x 12"
Net Weight	25 lbs
Material on product	Aluminum

Warranty

Journey Health & Lifestyle warrants the SoftSecure Rotating Transfer Bench frame to be free from defects in materials, workmanship and assembly for twelve (12) months from the original purchase date. The warranty does not extend to non-durable components such as the rubber tips.



 **journey**
HEALTH & LIFESTYLE®

 **journey**
SoftSecure

ROTATING TRANSFER BENCH

If you have any questions, call us at our toll free number:

1-800-958-8324

or visit www.shopjourney.com