


TRAVELITE TRANSPORT CHAIR (TC005GY)



- | | |
|---------------------------------------|--|
| 1. Handle Bar Left- TC00501L | 11. Footrest Assembly (pr)- TC005012 |
| Handle Bar Right- TC00501 | 12. Rear Wheel Axle & Spacer- TC005014 |
| 2. Handgrip- TC00502 | 13. Caster Axle- TC005015 |
| 3. Upholstery- TC00503 | 14. Cross Frame- TC005017 |
| 4. Arm Assembly Left- TC00504L | 15. Front Frame Tube Right- TC005018 |
| Arm Assembly Right- TC00504R | 16. Front Frame Tube Left- TC005019 |
| Armrest Tube- TC005013L | 17. Fork- TC005020 |
| 4a. Armrest Latch- TC00507 | Fork Bearing- STDS1060 |
| 4b. Arm Pad 10" Pad- TC00508 | 18. Rear Wheel- TC005021R |
| Arm Pad 8" Pad- TC005028 | 18a. Rear Wheel Bearing (pr)- STDS4S3602 |
| 5. Folding Hinge Left- TC00505L | 19. Caster- TC005021F |
| Folding Hinge Right- TC00505 | 19a. Caster Bearing- STDS1011 |
| 6. Telescoping Tube- TC00506 | 20. Wheel Lock Left- TC005022 |
| 7. Cross Frame Block (4 ea)- TC005011 | Wheel Lock Right- TC005023 |
| 8. Upholstery Strap- TC005027 | 21. Tipping Sleeve- STDS3J2422 |
| 9. Footrest Connector- TC00509 | 22. Frame Plug- TC005026 |
| 10. Footrest Tube- TC005010 | 23. Carry Bag- TC005024 |